Perkins Mk2

Hand-Held
Applanation Tonometer

Instruction Manual

Manuel d'instructions

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Manual de instrucciones

Please read these instructions
before using the equipment

CE 0120

	Page
English	1
French	34
German	67
Italian	101
Spanish	134

Contents

1. Introduction	2
1.1 Technical Description	2
1.2 Purpose of Use	3
1.3 Brief Description of the Appliance	3
2. General Layout	4/5
3. Safety and Regulatory Information	6
3.1 Safety Notices	6
3.2 Measuring Prism Replacement	7
3.3 Transportation and Storage	8
3.4 Notes on Usage	8
3.5 Statutory Requirements	9
3.6 Caution	9
3.7 Regulatory Information: CE Mark	9
3.8 Environmental Conditions of Use	10
3.9 EMC (Electromagnetic Compatibility) Declaration	11
3.10 Product Labelling and Safety Symbols	12
3.11 Safety Symbols	12
4. Taking a Pressure Measurement	13/14
5. To Prepare the Prism for Use	15
5.1 Remove Perkins Mk2 Tonometer From Case	15
5.2 Doubling Prism	15
5.3 To Insert Prism into Prism Holder	15
5.4 Tonosafe Disposable Prism	15
5.5 Batteries	15/16
6. Operating the Equipment	17
6.1 Preparing the Patient	17
6.2 Instruction to the Patient	18
7. Cleaning the Equipment	19
8. General Notes	20
8.1 Astigmatism	20/21
8.2 To Use the Instrument on Animals	21
8.3 Special Precautions	21
9. Measurement	22/23
9.1 Sources of Error	24
9.2 Inaccurate Horizontal Centering	25
9.3 Inaccurate Vertical Centering	26
9.4 Inaccurate Pressure	27
10. Case Contents	28
11. Calibration Check	29
12. Storage of Prism and Tonometer	30
13. Technical Specifications	31
References	32
Accessories	33

Whilst the information is provided in good faith and is based on the latest information available at the time of issue, this booklet gives only a general indication of product capacity, performance and suitability. Such information must not be taken as establishing any contractual or other commitment on the part of the manufacturer and in no way should be construed as a warranty or representation concerning the product.

1. Introduction

We would like to thank you for your decision to purchase this CLEMENT CLARKE product.

If the instructions in this manual are carefully followed, we are confident that this product will give you reliable and trouble-free usage.

We will be forced to disclaim any warranty and liability if the instrument is altered in any way or if routine maintenance is neglected or not carried out according to factory specifications.

Only original CLEMENT CLARKE accessories are to be used to ensure continued reliability of the instrument.

Because the **Perkins** tonometer is capable of being used in all positions, the mass of the Goldmann doubling prism must be controlled accurately. Only **HAAG-STREIT prisms should be used** or the disposable **Tonosafe** prism can be used to prevent cross-infection. Any other prisms may give a false reading if used in the **Perkins** tonometer.

1. 1 Technical Description

A pair of spiral wound flat springs connected in series are progressively tensioned by rotation of the operator's calibrated thumb-wheel. The resulting force stored in the spring is transmitted through a balanced arm to the applanation doubling prism or the Tonosafe prism mounted at its end.

The battery charger conforms to Class II insulation and has also been tested to EN60601.

This equipment is designed and manufactured to be as safe as is reasonably practicable when installed and used correctly by persons qualified to do so in accordance with the manufacturer's instructions. Any misuse, incorrect assembly or installation could result in the equipment becoming unsafe. It must not be repaired except by competent persons under the manufacturer's instructions.

1. Introduction *- continued*

1. 2 Purpose of Use

Applanation Tonometers are used at room temperature in the examination, diagnosis and documentation of the human eye.

Applanation tonometry may only be performed by healthcare professionals who are qualified by their respective authorities to do so.

1. 3 Brief Description of the Appliance

The Applanation Tonometer functions according to the “Goldmann” method: The measuring of the pressure required to maintain a uniform applanation of the surface of the cornea.

Advantages

- High accuracy of the measurement - the average deviation in any one single examination does not exceed ± 0.5 mm Hg.
- The intraocular pressure in mm Hg is found by multiplying the drum reading by ten.
- Scleral-rigidity is not to be taken into account as the small volumetric displacement of 0.56 mm^3 increases the intraocular tension by about 2.5% only.
- Massage effect DOES occur if repeated measurement occurs.
- There are no standardization and calibration difficulties.

Hospital and Optician use only

2. General Layout

Fig 1


- 1 Forehead rest
- 2 Forehead rest locking screw
- 3 Forehead rest stem
- 4 Prism holder
- 5 Filter Holder Assemblies
- 6 Illumination bulbs

2. General Layout - continued

Fig 2

- 7 Viewing lens
- 8 Milled thumb-wheel
- 9 Scale
- 10 Battery handle securing buttons
- 11 Battery handle

3. Safety and Regulatory Information

3. 1 Safety Notices

No user serviceable parts inside.

Do not open the Tonometer case moulding under any circumstances.

 **Caution:** In the event of any malfunction, contact your local distributor or visit our website to contact us directly.

This product is intended for indoor use only.

This product is not suitable for operation in the presence of a flammable anaesthetic mixture with air or with oxygen or with nitrous oxide.

Not protected against the ingress of water or other fluids.

The following instructions should be observed to ensure the safety of patients:

Avoid examinations in cases of eye infections or injured corneas

- There is no known contraindication for performing tonometry. Medical assessment is necessary and due care is to be exercised in cases of eye infections and injuries to the cornea.

Only clean and disinfected measuring prisms should be used

- Examinations should be conducted with cleaned and disinfected measuring prisms.
- Tonometer prisms are not dispatched in a disinfected state and should always be disinfected before use.
- The disinfecting of contact glasses e.g. with cotton wool or a damp cloth without putting them into disinfecting fluid is insufficient and should not be allowed.
- Incorrect medical disinfecting can cause patients to be infected and can also cause damage to the eyes.
- Residue from the disinfecting may cause a caustic reaction to the patient's eye, therefore rinse thoroughly.

3. Safety and Regulatory Information - *continued*

3. 1 Safety Notices - *continued*

No damaged prisms should be used

Before each use the contact surface of all prisms should be checked for contamination or damage (such as scratching, rips or sharp edges). A slit lamp microscope or Loupe is best used for this purpose, at a magnification of 10x to 16x.

- If cracks form, the disinfectant may get into the hollow interior and can enter the patient's eye while tonometry is taking place, leading to irritation of the cornea. Defective measuring prisms should not be used, to avoid any damage to the patient's cornea.

3. 2 Measuring Prism Replacement

- Due to the number of variables involved - type and concentration of disinfecting material used, number of patients, handling etc. - it is almost impossible to put a figure on the number of times and/or period of time that a measuring prism can safely be used for.
- Our experience, and that of our users, indicates that under 'normal' usage (i.e. as per the instructions contained within this instruction manual) a maximum of two years provides an adequate safety factor.
- HAAG-STREIT therefore recommend strongly that no measuring prism is used beyond two years from it's date of first usage (a measuring prism may be retained as a spare indefinitely - deterioration does not start until it is first used).
- The recommended time is overruled if the prism is in any way damaged within this time period - **damaged prisms should be immediately discarded.**

3. Safety and Regulatory Information - *continued*

3.3 Transportation and Storage

- The Applanation Tonometer conforms to ISO 15004-1:2006 for Storage and Transport Conditions and should always be stored in its carry case when not in use.

Transport Conditions

Temperature:	-40°C to +70°C
Relative Humidity:	10% to 95%
Atmospheric Pressure:	500hPa to 1060hPa
Vibration, Sinusoidal:	10Hz to 500Hz: 0.5g
Shock:	30g, duration 6 ms
Bump:	10g, duration 6 ms

Storage Conditions

Temperature:	-10°C to +55°C
Relative Humidity:	10% to 95%
Atmospheric Pressure:	700Pha to 1060hPa

- The Applanation Tonometer should be checked for damage after unpacking.
- Defective equipment should always be returned in appropriate packaging.
- The equipment may only be operated with the listed accessories (see page 33) according to the corresponding operating instructions.

3.4 Notes on Usage

- Only qualified and trained personnel may operate the equipment.
- The training of the operating personnel is the responsibility of the equipment user.
- This type of sensitive measuring equipment should be checked after exposure to external force (e.g. unintentional shocks or being allowed to fall) and if necessary, should be returned to the manufacturer for repair.
- The manufacturers of the equipment will not be liable for any loss or damage resulting from unauthorized intervention.

In particular, it is forbidden to loosen any screws; this may reduce the accuracy of the measurement. All ensuing warranty claims will be deemed null and void as a result of unauthorized intervention.

3. Safety and Regulatory Information - *continued*

3. 5 Statutory Requirements

- The Applanation Tonometer is designed as a Class IIa device under the Medical Device Directive 93/42/EEC.
- The Applanation Tonometer complies with the Tonometer regulation ISO 8612:2009 & ISO 15004-1: 2006.

3. 6 Caution

- The safety regulations displayed in the operating instructions are to be observed with special care.

Caution

Strictly observe all warning notices!

3. 7 Regulatory Information: CE-Mark

The products fulfil the requirements according to the following standards:

Safety: EN60601-1 Part 1: General Requirements for Safety, Medical Electrical Equipment

Electrical: EN60601-1-2 Part 2: Collateral Standard: Electromagnetic Compatibility

3. Safety and Regulatory Information - *continued*

3. 8 Environmental Conditions of Use

Temperature:	+10°C to +35°C.
Relative Humidity:	30% to 90%
Atmospheric Pressure:	800hPa to 1060hPa
Shock (without packing):	10g, duration 6 ms

3. Safety and Regulatory Information - *continued*

3. 9 EMC (Electro Magnetic Compatibility) Declaration

With the increased number of electronic devices such as PC's and mobiles, medical devices in use may be susceptible to electromagnetic interference from these devices. Which may result in incorrect operation of the medical device and create a potentially unsafe situation. Medical devices should also not interfere with other devices.

In order to regulate the requirements for EMC with the aim to prevent unsafe product situations, the EN60601-1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices.

Medical devices manufactured by Haag-Streit-UK Ltd conform to this EN60601-1-2 standard for both immunity and emissions.

Nevertheless, special precautions need to be observed:

The use of accessories and cables other than those specified by Haag-Streit-UK Ltd, with the exception of cables sold by Haag-Streit-UK Ltd as replacement parts for internal components, may result in increased emission or decreased immunity of the device.

The medical devices should not be used adjacent to or stacked with other equipment.

If adjacent or stacked use is unavoidable, the medical device should be seen to operate normally as it should in this situation.

Further guidance regarding the EMC environment (in accordance with EN60601-1-2) in which the device should be used is available at:

<http://www.haagstreituk.com/medicalquality>

3. Safety and Regulatory Information - *continued*

3.10 Product Labelling and Safety Symbols

The labels affixed to the equipment are there to remind you of safety precautions. Do not remove these labels.

3.11 Safety Symbols

The Perkins Mk2 Tonometer and its labels include the following markings:

	Caution (Consult 'Instructions for Use')
	This product complies with the essential requirements of the medical devices directive (93/42/EEC). Compliance has been verified by UK notified body per 0120 (SGS United Kingdom Ltd).
	Do not dispose of with household refuse. Please contact your local authority for your nearest recycling centre.
	Manufacturer
	Medical Electrical Equipment Type B
	Year of Manufacture
e.g. SN0001 T3	Serial number Perkins Mk2 Tonometer (Example)
	Caution of Hot surfaces - Do not touch

4. Taking a Pressure Measurement

The cornea is applanated with a measuring prism of organic glass that is fitted in a ring shaped holder at the end of the feeler arm. The contact surface of the measuring prism has a diameter of 7.0mm. It is planar with rounded edges to obviate any damage to the cornea.

The measuring prism is brought into contact with the eye. The measuring drum is turned to increase pressure on the eye until a continuously regular surface of 3.06 mm $\varnothing = 7.354 \text{ mm}^2$ has been applanated.

Position of the measuring drum	Force	Pressure
1	9.81 mN	1.33 kPa
		=10mm Hg

The relationship between the measuring drum position, force and pressure on the applanation surface

The intraocular pressure in mm Hg is found by multiplying the drum reading by ten.

SI-Units according to standard ISO 1000

According to ISO 1000 the derived unit of force is the N (Newton) and its decimal sub-multiple the mN (milli Newton).

The conversion of mm Hg into the SI pressure unit kPa (Kilo Pascal) can be made on the comparison scale in Fig 3.

The measurement of the flattened surface is made directly on the cornea. The built-in duplication system in the measuring prism divides the image and displaces the two semi circular halves from each other by 3.06 mm.

Before measuring, it will be necessary to anaesthetise the cornea locally, to place a strip of fluorescein paper in the conjunctival sac and switch on the blue filter of the tonometer.

The inner border of the ring represents the line of demarcation between the cornea flattened by applanation and the cornea not flattened.


Fig 3
Comparative scale

4. Taking a Pressure Measurement - *continued*

The major advantage of applanation tonometry is the small bulbous deformation that amounts to only 0.56 mm^3 . The values found by this method of tonometry are only slightly influenced by scleral rigidity and radius of corneal curvature.

The principle of applanation tonometry is simple. The careful construction of the apparatus is a guarantee of its continued correct performance.

Nevertheless, in order to obtain exact results, the method of use indicated must be strictly followed.

The readings obtained are based on the assumption that a cornea with a 'normal' thickness is being measured - alteration of the corneal thickness leads to changes in the measured IOP.

A 'normal' corneal thickness is considered to be within the area of 530-560 microns. There is no universally agreed formula for the adjustment of the measured tonometry reading against altered corneal thickness.

The users should bear these two facts in mind and keep themselves informed, through the current literature, with regard to the various formulae for adjusting measured tonometry readings in the event of thinned, or thickened, corneas being measured.

5. To Prepare the Prism for Use

5. 1 Remove Perkins Mk2 Tonometer from case

Remove the **Perkins** tonometer from the carrying case by holding the main body. Under no circumstances should it be lifted by the prism holder (**4** Fig 1).

5. 2 Doubling Prism

The doubling prism should be prepared for use by cleaning following the cleaning instructions supplied with the prisms.

5. 3 To Insert Prism into Prism Holder

Fig 4

To insert the prism into the prism holder (**4** Fig 1) gently press prism cone downward into the prism holder (**A** Fig 4) and push into position along horizontal axis (**B** Fig 4).

N.B. Failure to follow these instructions could result in severe damage being caused to the instrument.

5. 4 Tonosafe Disposable Prism

Alternatively, Tonosafe disposable prisms may be used. These are supplied as sterile **single use** items to help eliminate the risk of cross infection.

Tonosafe and the prism holder have the same mass as the original doubling prism, so no conversion of results is required. The disposable prism can be inserted in the same manner as above.

5. 5 Batteries

The battery handle will accept four standard size AA batteries. To fit them in the instrument proceed as follows:-

Press the two battery handle securing buttons (**10** Fig 2) simultaneously and withdraw the battery case (**11** Fig 2), which forms the handle portion of the tonometer. Place the batteries into the battery handle (**11** Fig 2) ensuring they are positioned as indicated by the symbols on the base (Fig 5a, p16).

5. To Prepare the Prism for Use - *continued*

5. 5 Batteries - *continued*

Replace the handle and advance the thumb-wheel to check that both bulbs are switched on.

Fig 5a

Fig 5b

Note:

Ensure batteries are inserted into handle with the correct polarity, as marked on the base and in accordance with the Instruction Manual (Page 16, Batteries). Failure to do this will cause short circuiting of batteries and dangerous overheating will occur.

Disposal of used batteries

Discharged non-rechargeable batteries should not be incinerated, but should be disposed of in accordance with local environmental legislation.

Note: *All batteries are subject to slow deterioration. Under no circumstances should batteries be left indefinitely in the instrument as this could result in the instrument becoming damaged. They should be removed after use and stored in the carrying case.*

Bulbs

The bulbs, two in number, are “lens” lamps (**6** Fig 1). They are clear and require the blue filter sleeves to be fitted over them (**5** Fig 1). The bulbs are fixed in position to create maximum fluorescent effect during measurement. Bulbs must be screwed fully and firmly in and periodically inspected for discolouration.

Forehead Rest

This is the concave disc (**1** Fig 1) which snaps on to the rounded end of the forehead rest stem (**3** Fig 1). To adjust the stem, loosen the locking screw (**2** Fig 1), slide the stem to the required position as described in Method of Use and tighten the locking screw. The adjustment allows for considerable variation in the size and shape of a subject’s physical features. The disc of the forehead rest should be cleaned regularly.

6. Operating the Equipment

Caution

It is imperative to read chapter 3. Safety and Regulatory Information, and to observe its precautions before operating the equipment.

Perfect functioning can only be guaranteed with the use of an original HAAG-STREIT prism or Tonosafe.

Sterilized fluorescein paper strips should always be used, as pathogenic excitors thrive well in fluorescein solutions.

Regular calibration of the Tonometer is advised.

Do not use without prism fitted.

6. 1 Preparing the Patient

1. Both eyes must always be anaesthetized (e.g. 2-3 drops each of an anaesthetic within half a minute) to reduce movements of the lids during examination.
2. Place a fluorescein paper strip near the external canthus in the lower conjunctival sac. After a few seconds the lacrimal fluid is sufficiently coloured and the paper can be removed.

Alternatively combined anaesthetic fluorescein eye drops can be used.

When using drops, a solution of sodium fluorescein 0.25% to 0.5% is advisable.

Caution: If the prism touches the cornea without any force applied, the pressure arm can vibrate and may disturb the patient. For this reason, the scale is set to read 1 before applying the prism face to the cornea.

The **Perkins** tonometer should be held so that the thumb rests on the milled wheel (Fig 2, **8**) controlling the spring.

If the forehead rest (Fig 1, **1**) is to be used, the stem (Fig 1, **3**), should be extended after loosening the locking screw (Fig 1, **2**). When the forehead rest is correctly positioned, tighten the locking screw and the instrument can then be brought towards the patient so that the face of the prism cone comes into contact with the cornea.

It is usually easier to hold the tonometer obliquely with the handle slanted away from the nose (Fig 6). Care should be taken to prevent the prism touching the eyelid margin.

6. Operating the Equipment - *continued*

6. 1 Preparing the Patient - *continued*

When seen through the viewing lens (Fig 2, **7**) the semicircles of fluorescein should now be clearly visible and appear equal in size, denoting that the contact point is central on the cornea.

The light is switched off by turning the thumb-wheel until the scale reading (Fig 2, **9**) is below zero.

! **Note:** *The lighting is designed to give optimum fluorescence only when the prism cone is at its calibrated position i.e. when the rim of the prism holder is level with, or not more than 1 mm above the case surface.*

6. 2 Instruction to the Patient

Fig 6

Instruct the patient to look straight ahead or slightly upward and if necessary use a fixation light for this purpose.

It is recommended that patients are repeatedly asked to keep their eyes wide open during the examination. If required, the lids of the examined eye may be held open by the examiner's fingers provided that no pressure is applied to the eye.

7. Cleaning the Equipment

Cleaning (manual) – non-immersion for Ophthalmic medical devices.

Scope of Action

Non-immersion manual cleaning methods are appropriate for low risk items (those items that come into contact with intact skin or do not contact the patient) where soaking in aqueous solution, e.g. electrical and electronic equipment, will compromise the device.

Alcohol wipes should be used to clean electrical contacts on equipment.

Equipment required:

1. A warm water/detergent solution at correct/recommended dilution.
2. A clean, disposable, absorbent, non-shedding cloth for application of detergent solution.
3. A clean, disposable, absorbent, non-shedding cloth for drying equipment.
4. An appropriate chemical neutraliser, first aid kit and eyewash, in case of splashing with detergent.

Procedure:

Ensure that it is disconnected from the charger and the unit is switched off (pointer reaches just below zero) before commencing the cleaning procedure.

Wearing protective clothing, immerse the cleaning cloth in the detergent solution and wring thoroughly.

Commencing with the upper surface of the item, wipe thoroughly ensuring that detergent solution does not enter the electrical components.

Periodically rinse the cloth in clean water and repeat the above steps.

Surfaces should be carefully hand-dried using a fresh dry cloth.

Note: Non-immersion, manual cleaning is not a disinfection process, but where an alcohol wipe is used to dry surfaces, this may have a disinfecting effect.

Safely dispose of cleaning materials and alcohol wipes, if used.

The above was written with reference to: Sterilization Disinfection and cleaning of Medical Devices & Equipment: Guidance on Decontamination from the Microbiology Advisory Committee to the Department of Health Medical Devices Agency (MAC manual section 2 page 20 updated April 2005)

8. General Notes

Note

Only the shortest possible measuring time periods should be carried out on each eye. Should drying of the cornea epithelium be observed, then vision and field of sight should be checked in advance.

The measuring procedure can be repeated several times. Excited or anxious patients often display a higher ocular pressure during the first measurement.

For this reason, a reduction of tension is experienced during the first few minutes, as patients notice that tonometric examinations cause no unpleasant effects.

A test measuring procedure should therefore be first made on each eye. These results can be disregarded. Thereafter, three measuring procedures are to be carried out on each eye. The readings will be correct when the pressure has been stabilised. When measuring procedures have been correctly carried out, divergences in the results will be only of the order of ± 0.5 mm Hg.

When the measuring procedure lasts too long on one eye, drying will occur on the cornea epithelium of both eyes. A ring of fluorescein-positive deposits around the contact positions of the cornea and the measuring element will occur on the eye under examination.

On the other eye, map-like fluorescein-positive dried places will occur, which will not allow a reliable measuring procedure.

Extensive drying phenomena disappear after a short time without treatment becoming necessary. Focused vision will be affected by such fine epithelium defects.

8. 1 Astigmatism

If the cornea is spherical, measurements can be made on any meridian, but it is most convenient to do it on the 0° meridian.

This is not so when eyes with higher corneal astigmatism that 3 dioptres are examined, as the flattened areas are not circular but elliptic.

It has been calculated that, in cases of larger corneal astigmatisms, a surface of 7.354 mm^2 ($\varnothing 3.06 \text{ mm}$) is to be applanated, when the measuring prism is at an angle of 43° to the meridian of the greatest radius.

8. General Notes - *continued*

8. 1 Astigmatism - *continued*

For example

If the corneal astigmatism amounts to

$$\begin{aligned} 6.5 \text{ mm}/30^\circ &= 52.0 \text{ D}/30^\circ \text{ and} \\ 8.5 \text{ mm}/120^\circ &= 40.0 \text{ D}/120^\circ \end{aligned}$$

the graduation value 120° of the prism is set at the 43° mark of the prism holder.

If there is a corneal astigmatism of

$$\begin{aligned} 8.5 \text{ mm}/30^\circ &= 40.0 \text{ D}/30^\circ \text{ and} \\ 6.5 \text{ mm}/120^\circ &= 52.0 \text{ D}/120^\circ \end{aligned}$$

the graduation value 30° is set at the 43° mark.

In other words, set the axial position of the greatest radius, that is the axis of a minus cylinder, on the prism graduation at the red mark on the prism holder.

8. 2 To use the instrument on animals

Unfortunately, most animal eyes have a cornea that behaves differently from the human cornea and the Fick-Maklakov law does not apply for the standard diameter of applanation. It is necessary to remove the patient's forehead rest and calibrate the tonometer for each species used. This can be done quite simply by connecting the anterior chamber to a reservoir and taking tonometer readings at different pressure levels. In most cases, it will be found that the tonometric reading is lower than that obtained for the human eye. The eyes of rhesus and cynomolgus monkeys behave very similarly to the human eye, but those of the rabbit and cat give lower readings.

8. 3 Special Precautions

The **Perkins** hand-held tonometer is a delicate precision instrument, carefully inspected and certified. Each instrument has been tested at 0 - 5 gram positions on the scale (1 gram increments) and found to be correct within the limits specified on the certificate accompanying each tonometer. It is desirable to clean the calibrating weights, as well as the prism, in an antiseptic solution, before checking the instrument, to minimize the risks of cross-infection (see cleaning instructions accompanying the prisms for decontamination instructions).

9. Measurement

Measurement of the force necessary to flatten a corneal surface of constant (\varnothing 3.06mm) size.


Fig 7

Adjust the force by turning the thumb-wheel (Fig 2, **8**) until the inner margins of the semicircles coincide (Fig 7).

The tonometer is removed from the eye and the reading noted. The large divisions of the scale represent grams and the small divisions 0.2 gram.

The reading is multiplied by ten to give the tension in millimetres of mercury (mm Hg). Readings should be repeated until a steady value is obtained; this will usually be slightly lower than the initial reading.

If the semicircles appear large and are not reduced by diminishing the force of the spring, the tonometer has been pushed too close to the eye and withdrawing it slightly will bring the prism within the range of free movement.

If it is desired to use the **Perkins** tonometer without any forehead rest, the instrument can be steadied by resting the fingers against the patient's malar region while the prism is in contact with the cornea.

9. Measurement - continued

1. Immediately before taking the measurements, the patient should be made to close the eyes briefly so that the cornea becomes sufficiently moistened with the lacrimal fluid and fluorescein.

2. By moving the tonometer forward, the measuring prism is brought into contact with the centre of the cornea over the area above the pupils. The limbus of the cornea will be illuminated with a bluish hue. This illumination is best observed by direct sight from the opposite side of the illumination unit by the examiner.

As soon as the limbus of the cornea illuminates, cease any further forward movement of the tonometer immediately.

3. After obtaining contact, commence observing the cornea. A regular pulsation of the two semi circular fluorescein rings, which can be of varying size dependent on the ocular pressure, when the measuring drum is at setting 1, will show that the tonometer is in the correct measuring position.

4. The pressure on the eye is increased by turning the milled thumb wheel on the tonometer, until the edges of both fluorescein rings just meet (Fig. 9). The edges cross over each other with the pulsation of the eye. The width of the fluorescein ring around the contact position of the measuring prism should be $1/10$ approx. of the diameter of the applanation surface (0.3 mm).

5. Reading the scale:

- The reading
- multiplied by 10
- gives the ocular pressure in mm Hg

Note

The following are schematically simplified illustrations!


Fig. 8

Semi circular images in the centre of the field of sight


Fig. 9

Correct final position

9. 1 Sources of Error

4.5.1 Fluorescein ring / Distance to the patient





Result	Image in ocular	Cause	Correction
The fluorescein ring is too wide		The measuring prism was not dried after cleaning, or the eyelids came into contact with the measuring prism whilst measuring.	The tonometer should be withdrawn and the measuring prism dried with a cotton wool swab.
The fluorescein ring is too narrow		The lacrimal fluid has dried.	Let the patient close the eyes once or twice. Then repeat the measuring procedure.
No semi circular rings appear, only the centre line		The measuring prism is not touching the cornea. Should the patient withdraw their head slightly, the irregular pulsations will occur and the prism will only contact the eye intermittently. Should the patient withdraw the head further, then the fluorescein rings will disappear altogether.	Steady patient's head.
Both of the too large semi circular rings appear partly		Should the tonometer be moved forward towards the patient, or should the patient move towards the tonometer whilst the measurement is being taken, then the feeler arm will come into contact with a sprung stop piece. The applanation surface is then too large.	The image will not change when turning the milled thumb wheel. Withdraw the tonometer until the regular pulsations of a corresponding smaller applanation surface indicate the correct measuring position and pressure changes lead to immediate applanation surface changes.

Fig 10

9. 2 Inaccurate Horizontal Centering






Result	Image in ocular	Cause	Correction
The upper semi circular ring appears partly		The measuring prism is not focused on the eye. The eye is too far on the right.	Move the tonometer to the right.
The upper semi circular ring appears completely - the lower ring partly		The measuring prism is not focused on the eye. The eye is still too far on the right.	Move the tonometer to the right.
Correct position: The edges of both semi circular rings meet exactly in the centre of the ocular		The measuring prism is focused precisely on the eye.	
The lower semi circular ring appears completely - the upper ring partly		The measuring prism is not focused on the eye. The eye is still too far on the left.	Move the tonometer to the left.
The lower semi circular ring appears partly		The measuring prism is not focused on the eye. The eye is too far on the left.	Move the tonometer to the left.

Fig 11

9. 3 Inaccurate Vertical Centering





Result	Image in ocular	Cause	Correction
A semi circular ring appears partly in the upper half		The measuring prism is not focused on the eye. The eye is too far up.	Move the tonometer upwards.
The ring appears completely in the upper half		The measuring prism is not focused on the eye. The eye is still too far up.	Move the tonometer upwards.
The ring appears almost completely in the upper half, and partly cut in the lower half		The measuring prism is not focused on the eye. The eye is still too far up.	Move the tonometer upwards.
Two partly cut rings appear, the larger one in the upper half		The measuring prism is not focused on the eye. The eye is still too far up.	Move the tonometer upwards.
Correct position: The edges of both semi circular rings meet exactly in the centre of the ocular		The measuring prism is focused precisely on the eye.	


Fig 12

9. 4 Inaccurate Pressure

Result

The outer edges of the fluorescein rings touch each other

Image in ocular


Cause

Pressure is too strongly reduced

Correction

Increase the pressure by turning the milled thumb wheel.


The fluorescein rings coincide and form a line

The pressure is reduced

Increase the pressure by turning the milled thumb wheel.


The fluorescein rings do not touch each other

The pressure is too strongly increased

Reduce the pressure by turning the milled thumb wheel.

Correct position: The inner edges of the fluorescein rings touch each other

Accurate pressure - precise focusing of the measuring prism

Fig 13

10. Case Contents

- | | | | |
|---|--------------------------------------|---|------------------|
| 1 | Perkins Tonometer Mk2 | 5 | Five-gram weight |
| 2 | Forehead Rest | 6 | Batteries |
| 3 | Doubling Prism/Black Tonosafe Holder | 7 | Spare Bulb |
| 4 | Two-gram weight | 8 | Setting Block |

Fig 14

11. Calibration Check (Fig 15)

Regular calibration of the Tonometer is advised.

 **Please see 7.3 Special Precautions**

The calibration of the tonometer should be checked at regular intervals or at any time when the practitioner feels that an immediate spot check is desirable.

Proceed as follows:-

1. Remove the battery handle (**11** Fig 2) from the tonometer by pressing the two handle securing buttons (**10** Fig 2) and withdrawing it from the body portion. Insert the prism into the holder (**4** Fig 1).

Check for zero position

2. Adjust the milled thumb-wheel (**8** Fig 2) so that the scale reading is below zero by the full thickness of one scale line. Whatever the attitude of the instrument in the hand the prism should tend towards its backward position.

3. Set the scale reading above zero by the full thickness of one scale line. The prism should tend towards its forward position regardless of the attitude of the tonometer.

Check at the 2 g position

4. To make a rough check adjust the milled thumb-wheel (**8** Fig 2) so that the scale reading is below the 2 mark by the full thickness of one scale line. Place the instrument on a flat horizontal surface with the setting block under the body and the prism cone upwards as shown in Fig 16. Place the 2 g weight, recess downward, centrally on the prism cone and check that it carries the prism positively down to its lowest position.

5. Also, as a rough check, set the scale reading above the 2 mark by the full thickness of the scale line and place the instrument on the setting block as before. Place the 2 g weight on the prism and check that it remains in its highest position.

(If the prism is accidentally pushed down to its lowest position with the weight on the prism, it should rise again to its highest position and freely remain there).


Fig 15

If desired, the **Actual** point of balance can be established by adjusting the milled thumb-wheel (**8** Fig 2) until the prism floats at the calibration-position (*rim of prism holder level with, or not more than 1 mm above the case surface*).

11. Calibration Check - continued

Check at the 5 g position

6. Repeat the process as described in 4. and 5. above, but using the 5 mark on the scale and a setting one-and-a-half lines width above and then below the 5 mark and with the 5 g weight on the prism cone. The actual point of balance can again be found if required.

Note: *The tolerances in the zero, 2 g and 5 g tests above represent ± 0.05 , ± 0.05 and ± 0.075 g respectively. If the tonometer is found to be outside the stated limits it should be returned to your supplier for attention.*

12. Storage of Prism and Tonometer

The doubling prism must always be removed from the tonometer, cleaned following the cleaning instructions supplied with the prisms and inserted, cone downward, into the hole in the foam interior of the carrying case (Fig 14, **3**). The tonometer should then be returned to the carrying case.

Do not pack the **Perkins** tonometer with the prism still in the instrument.

Never store the instrument with the spring under load - i.e. return the milled thumb-wheel to below zero.

The suppliers cannot accept responsibility for the performance, reliability or safety of this equipment if it is used, serviced or modified by unauthorized persons.

Exterior cleaning of the tonometer body should be carried out using an antiseptic wipe, a cloth moistened in a mild soap solution or with 70% isopropyl alcohol.

Solutions containing alcohol, benzene etc., should not be used.

Note

Repairs

In the event of damage do not attempt local repairs but return the tonometer in its case, carefully packed, to your manufacturer for attention.

Warranty

Unauthorized attempts at repair will immediately nullify the warranty.

13. Technical Specifications

Measuring range

0 - 50 mm Hg

Measurement uncertainty of the force impingement on the measuring prism over a measuring range from 0 - 49.03mN

Standard divergence:
 $0.49\text{mN} \leq 3s \leq 1.5\%$ of rated value

Net Weight, packed

1.0kg

References

- Perkins ES:** Hand-held applanation tonometer.
Br J Ophthalmol 1965 **49** 591.
- Perkins ES:** Tonometry in experimental animals.
Evaluation of Drug Effects on the Eye 1968 43-47.
- Heid S and Flament J:** Aplanotonométrie à main.
Bulletin des Sociétés d'Ophtalmologie de France 1968 **68** 655.
- Flament J:** Aplanotonométrie à main chez le sujet couché.
Bulletin des Sociétés d'Ophtalmologie de France 1969 **69** 640.
- Whitty HPB:** Trial results of a hand-held applanation apparatus.
Br J Ophthalmol 1969 **53** 664.
- Chignell AH:** Use of the Perkins hand-held tonometer in retinal detachment surgery.
Br J Ophthalmol 1971 **55** 644-646.
- Dunn JS and Brubaker RF:** Perkins applanation tonometer, clinical and laboratory evaluation.
Arch of Ophthalmol 1973 **89** 149-151.
- Radtke ND and Cohan BE:** Intraocular pressure measurement in the newborn.
Am J Ophthalmol 1974 **78** 501-504.
- Al-Abrak MH and Samuel JR:** Effects of general anaesthesia on the intraocular pressure in man.
Br J Ophthalmol 1974 **58** 806-810.
- Fujikawa LS, Salahuddin SZ, Palestine AG, Masur H, Nussenblatt RB and Gallo RC:** Isolation of human T-lymphotropic virus type III from the tears of a patient with the acquired immunodeficiency syndrome.
Lancet 1985 **2** 529-530.
- Recommendations for preventing possible transmission of human T-lymphotropic virus type III/lymphadenopathy- associated virus from tears.
Morbidity and Mortality Weekly Report 1985 **34** 533-534.

Catalogue Description

Number

5805001

Standard combination comprising - The Perkins hand-held tonometer, with extra bulb, doubling prism, patient's forehead rest, batteries, instruction manual and Certificate of Accuracy. The setting block and weights (5g and 2g) are also carried in the case provided.

Replacement components:

5805030

Carrying case

5801004

Five-gram weight

5801005

Two-gram weight

1146005

Bulb, 2.8V May lens

5805355

Filter Holder Assy

1158006

Alkaline battery, Type AA

5805009

Setting block

5805081

Patient's forehead rest

1902100

User's instruction manual

5801068

Doubling prism

Optional Accessories:

5805380

Rechargeable battery handle, complete

5805490

Battery charging unit

5814004

Perkins Examination Telescope Assembly

5805100

1 box of 100 prisms and 5 prism holders

5805080

1 box of 100 prisms and 5 prism holders. (USA only)

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Références

Perkins ES : Hand-held applanation tonometer. Br J Ophthalmol 1965 **49** 591.

Perkins ES : Tonometry in experimental animals.
Evaluation of Drug Effects on the Eye 1968 43-47.

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Chignell AH : Use of the Perkins hand-held tonometer in retinal detachment surgery.
Br J Ophthalmol 1971 **55** 644-646.

Dunn JS and Brubaker RF : Perkins applanation tonometer, clinical and laboratory evaluation. Arch of Ophthalmol 1973 **89** 149-151.

Radtke ND and Cohan BE : Intraocular pressure measurement in the newborn.
Am J Ophthalmol 1974 **78** 501-504.

Al-Abraak MH and Samuel JR : Effects of general anaesthesia on the intraocular pressure in man.
Br J Ophthalmol 1974 **58** 806-810.

Fujikawa LS, Salahuddin SZ, Palestine AG, Masur H, Nussenblatt RB and Gallo RC : Isolation of human T-lymphotropic virus type III from the tears of a patient with the acquired immunodeficiency syndrome.
Lancet 1985 **2** 529-530.

Recommendations for preventing possible transmission of human T-lymphotropic virus type III/lymphadenopathy- associated virus from tears.
Rapport hebdomadaire de morbidité et de mortalité 1985 **34** 533-534.

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Perkins ES: Hand-held applanation tonometer.

Br J Ophthalmol 1965 **49** 591.

Perkins ES: Tonometry in experimental animals.

Evaluation of Drug Effects on the Eye 1968 43-47.

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Arch of Ophthalmol 1973 **89** 149-151.

Radtke ND and Cohan BE: Intraocular pressure measurement in the newborn.

Am J Ophthalmol 1974 **78** 501-504.

Al-Abraak MH and Samuel JR: Effects of general anaesthesia on the intraocular pressure in man.

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Lancet 1985 **2** 529-530.

Recommendations for preventing possible transmission of human T-lymphotropic virus type III/lymphadenopathy- associated virus from tears.

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Whitty HPB: Trial results of a hand-held applanation apparatus.

Br J Ophthalmol 1969 **53** 664.

Chignell AH: Use of the Perkins hand-held tonometer in retinal detachment surgery.

Br J Ophthalmol 1971 **55** 644-646.

Dunn JS and Brubaker RF: Perkins applanation tonometer, clinical and laboratory evaluation.

Arch of Ophthalmol 1973 **89** 149-151.

Radtko ND and Cohan BE: Intraocular pressure measurement in the newborn.

Am J Ophthalmol 1974 **78** 501-504.

Al-Abrak MH and Samuel JR: Effects of general anaesthesia on the intraocular pressure in man.

Br J Ophthalmol 1974 **58** 806-810.

Fujikawa LS, Salahuddin SZ, Palestine AG, Masur H, Nussenblatt RB and Gallo

RC: Isolation of human T-lymphotropic virus type III from the tears of a patient with the acquired immunodeficiency syndrome.

Lancet 1985 **2** 529-530.

Recommendations for preventing possible transmission of human T-lymphotropic virus type III/lymphadenopathy-associated virus from tears.

Morbidity and Mortality Weekly Report 1985 **34** 533-534.

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Tel: +44 (0)1279 414969

Fax: +44 (0)1279 635232

Patented in Great Britain
& countries throughout the world

Instruction Manual no: 1902102
Issue no: 4 10/15

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Haag-Streit UK Ltd.,