INSTRUCTIONS FOR USE

ENGLISH

Purpose of use

The measuring prism is a reusable tonometer head, which serves as an accessory to the Haag-Streit tonometers AT 900, AT 900 D and AT 870.

1. Safety

1.1 Important general notes

WARNING!

- The measuring prisms listed in this instruction manual are not sterilized when shipped. They must be cleaned and disinfected prior to their first use.
- All the warnings and notes in the instruction manuals of the measuring prisms and Desinset must be followed without exception.
- Examinations should only be performed with undamaged, cleaned and disinfected measuring prisms.
- The measuring prisms used must not be damaged as they might otherwise damage the patient's cornea.
- Improper preparation can result in the transmission of diseases to the patient and user as well as damage to the measuring prisms.
- Residue from cleaning agents and disinfectants can irritate and burn the patient's eye.
- As a rule, the measuring prisms may be prepared together but should not be prepared with any other products.

1.2 Operation

WARNING!

- The operation and preparation must only be performed by qualified, trained personnel, whose training corresponds to the task of the operator.
- Appropriate professional assessment and caution are necessary.

1.3 Warranty and product liability

- The measuring prisms must be handled as described in the "Safety" section.
- Improper handling can damage the product. This voids all guarantee claims.
- Continuing to use a product damaged by incorrect handling can lead to personal injury. In this case, the manufacturer accepts no liability.

WARNING!

- Tonometer measuring prisms must be prepared in accordance with the specified process. This process is described in the separate instructions "Cleaning and Disinfection".
- Only use clean, residue-free and undamaged measuring prisms. Prior to every use, the contact surface of the measuring prism must be inspected for contaminants and damage (scratches, cracks and sharp edges). This must be performed with a slit lamp microscope at **10x to 16x magnification**.


2. Period of use

WARNING!

How long may measuring prisms be used?

Due to the large number of variables to be taken into account (type and concentration of the disinfectants used, number of patients, handling etc.) it is practically impossible to make a statement as to how often and/or how long a measuring prism can be safely used. The measuring prisms are stamped with an expiry date (□ YYYY-MM-DD). They must not be used beyond this date. Haag-Streit recommends a maximum service life of two years before the expiry date. This service life applies subject to normal conditions of use, i.e. in compliance with the instructions contained in these instructions for use. The period of use begins with the first application. The given time periods are not valid for damaged measuring prisms, they must be **immediately replaced**.

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HAAG-STREIT AG